

# COLLECTING SAMPLES FOR DNA TESTING IN SWINE




WASHINGTON STATE UNIVERSITY  
EXTENSION

There are three major sample collection methods for DNA testing in swine: tissue, blood, and hair bulbs. Ease and flow of sampling, storage of materials before and after collection, and cost are all important considerations for determining which method to use.



## TISSUE SAMPLING UNITS (TSU)

- Tissue sampling units, or TSUs, punch a small tissue sample from the ear into a tube with preservation solution.
- This sampling procedure can be combined with tagging, when the sample is taken at the same time as the tag is inserted into the ear.
- TSUs can be stored at room temperature in a cool, dry place before shipping at room temperature.





## WHOLE BLOOD & BLOOD CARDS

- Blood can be collected from a tail or jugular draw into a purple top tube (containing EDTA to prevent clotting).
- Whole blood should be stored in the fridge and shipped within a few days on ice packs.
- Blood can also be spotted onto a FTA blood card covering the entire circle or square on the card, and be sure to let the card dry before storing in a cool, dry place.
- Blood cards can be shipped at room temperature.



## HAIR BULBS

- Hair should be pulled and not clipped from the animal. First, comb out dirt and loose hairs from the area.
- Use a clean tool such as pliers to pull 20-30 hairs from the poll or neck area.
- Make sure to get the hair bulbs or roots, as these contain the DNA. Hair can be placed in hair cards with adhesive, stored in a cool, dry location and shipped at room temperature.



## SAMPLING TIPS

- Keep detailed and clear records during sampling to prevent duplicated or skipped animals.
- Make sure samples are not contaminated with manure, dirt, and debris.
- Clean sampling tools between animals to prevent cross-contamination of DNA.
- Store samples out of sunlight and heat, and ship accordingly to prevent DNA degradation.

**Kimberly M. Davenport**, Assistant Professor, Department of Animal Sciences, Washington State University  
**Sarah M. Smith**, Regional Animal Science Specialist, Washington State University Extension Moses Lake  
**Holly L. Neibergs**, Professor, Department of Animal Sciences, Washington State University

# REFERENCES AND ADDITIONAL RESOURCES

- DNA collection brings value to cattle producers. Scott Holt, Merck Animal Health. <https://www.merck-animal-health-usa.com/species/cattle/cattle-insights/beef/value-of-cattle-dna-collection>.
- DNA sample collection. Megan Rolf, Oklahoma State University and eBEEF.org. [https://beef-cattle.extension.org/wp-content/uploads/2019/09/2015-2\\_DNA\\_Sample\\_Collection\\_Arial.pdf](https://beef-cattle.extension.org/wp-content/uploads/2019/09/2015-2_DNA_Sample_Collection_Arial.pdf).
- DNA testing: News you need to know. Doug Newcom, National Swine Registry. [https://nationalswine.com/resources/docs/DNA\\_testing\\_story.pdf](https://nationalswine.com/resources/docs/DNA_testing_story.pdf).
- Pig DNA collection. UC Davis Veterinary Genetics Laboratory. <https://vgl.ucdavis.edu/sample-collection/pig>.
- SOP: Blood collection in swine. Virginia Tech University Veterinarian & Animal Resources. [https://ouv.vt.edu/content/dam/ouv\\_vt\\_edu/sops/large-animal/sop-swine-blood-collection.pdf](https://ouv.vt.edu/content/dam/ouv_vt_edu/sops/large-animal/sop-swine-blood-collection.pdf).
- Tissue sampling unit DNA collection. Neogen. <https://www.neogen.com/neocenter/blog/tsu-dna-collection/>.

For more information, visit: <https://wpnetwork.extension.wsu.edu/animal-genomics/>

## QUESTIONS? CONTACT US!

Dr. Kim Davenport

Email: [kimberly.davenport@wsu.edu](mailto:kimberly.davenport@wsu.edu)

Phone: (509) 335-4124

Website: <https://wpnetwork.extension.wsu.edu/animal-genomics>



WASHINGTON STATE UNIVERSITY  
EXTENSION

**Kimberly M. Davenport**, Assistant Professor, Department of Animal Sciences, Washington State University

**Sarah M. Smith**, Regional Animal Science Specialist, Washington State University Extension Moses Lake

**Holly L. Neibergs**, Professor, Department of Animal Sciences, Washington State University