

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL

USTUR 210: REMOVAL OF PLUTONIUM FROM SAMPLES TO BE ANALYZED FOR AMERICIUM ONLY

Purpose	Sample preparation for Am determination	Method Number	USTUR 210
Original Date	10/1/96	Author	Dorothy Stuit
Revision Number	0	Approved By	Jim Elliston
Revision Date	10/1/96	Approval Date	1/31/01

SAFETY NOTE: Before beginning this procedure, read all of the Material Safety Data Sheets for the chemicals listed in Section 3 of this procedure.

1. Principle of Method

- 1.1. This procedure is for removal of plutonium from samples being analyzed only for americium. This ensures that ^{238}Pu does not interfere with the ^{241}Am peak in alpha spectrometry measurements.
- 1.2. An aliquot of the tissue sample solution is selected for analysis based on wet sample mass (or bone ash weight) and estimated sample alpha activity.
- 1.3. Both Americium-243 and Plutonium-242 tracers are added to the aliquot. Americium is separated from plutonium by anion exchange.
- 1.4. The americium portion is purified by continuing with procedure USTUR 300: Anion Exchange Isolation of Americium from Prepared Tissue Solutions.
- 1.5. Chemical losses are identified and corrected on the basis of the recovery of the added ^{243}Am tracer.
- 1.6. Separation of Pu from Am is identified by absence of ^{242}Pu during alpha spectrometry.

2. Minimum Detectable Activity (MDA)

Not applicable.

3. Accuracy and Precision

Not applicable.

4. Apparatus

- 4.1. Hot plates.
- 4.2. Magnetic stirring hot plates.
- 4.3. Hot plate thermometer (100 - 150°C).
- 4.4. Watch glasses: assorted sizes to fit beakers used.
- 4.5. Bio-Rad Ion-exchange columns (Fig. 1): borosilicate glass barrel with polypropylene reservoir, column tip, stop cock, and bed support; 20 cm long by 1.0 cm i.d. or equivalent.
- 4.6. Rack: to support ion-exchange columns.
- 4.7. Glass beads: 2 or 3 mm diameter.
- 4.8. Graduated cylinders.
- 4.9. Beakers: various sizes.
- 4.10. Metric ruler.
- 4.11. Stir bars: Teflon-coated.
- 4.12. Wash bottles: 50 mL.

5. Reagents

- 5.1. Nitric acid (concentrated 69-71%, reagent-grade).
- 5.2. Nitric acid (8M). Add 500 mL concentrated HNO₃ to 400 mL of nanopure . Dilute to 1000 mL with nanopure water.
- 5.3. Sodium nitrite (reagent-grade).
- 5.4. Bio-Rad anion exchange resin (AG1-X4, 100-200 mesh) chloride form. Make a slurry of half resin, half nanopure water in a wash bottle.
- 5.5. Plutonium-242 tracer solution.
- 5.6. Americium-243 tracer solution.

6. Procedure

5.7. Sample preparation.

- 6.1.1 Refer to procedure USTUR 100, Tissue Ashing, Sample Dissolution, Sample Aliquot Selection, and Tracer Addition.

NOTE: Add both ^{243}Am and ^{242}Pu tracers (~10 dpm each).

6.2 Sample dissolution.

- 6.2.1. To the previously prepared sample, add 80 mL of 8 M HNO_3 . Use a stirring hot plate with a stir bar to ensure complete dissolution of the sample. Cover the beaker with a watchglass while heating (100-130°C).

6.3. Anion exchange separation.

- 6.3.1. Using a ruler, mark the ion exchange column 9 cm above the frit.
- 6.3.2. Fill the column with the slurry of AG1-X4 resin to a settled depth of 9 cm.
- 6.3.3. Rinse all of the resin down from the sides of the column using distilled water, then 8 M HNO_3 . Once all of the acid has drained, add approximately 1 cm depth of glass beads on top of the resin to prevent disturbance as reagents are added.
- 6.3.4. Wash the column with 45 mL of 8 M HNO_3 , collecting the solution in a waste beaker. Discard this solution into a hazardous waste car boy.
- 6.3.5. While conditioning the column, add approximately 0.5 g of sodium nitrite to the sample solution.
- 6.3.6. Swirl the beaker, cover with a watchglass, and heat the sample quickly, just to a boil, with a minimum amount of sample volume reduction.
- 6.3.7. Cool the sample to room temperature.
- 6.3.8. Place a 250 mL collection beaker under the column to catch the americium portion of the sample.
- 6.3.9. Pour the sample onto the column.
- 6.3.10. When the solution has drained to the top of the resin bed, rinse the empty sample beaker with 25 mL of 8 M HNO_3 , adding the rinse to the column.
- 6.3.11. Repeat step 6.3.10.

6.3.12. When the solution has reached the top of the resin bed, remove the americium collection beaker and replace it with a waste beaker. Evaporate the americium fraction to dryness on a hot plate (~120-150°C).

NOTE: Proceed to procedure USTUR 300: Anion Exchange Isolation of Americium from Prepared Tissue Solutions, for the Determination of Americium.

6.3.13. Rinse the column with 50 mL nanopure water. Discard the rinse into a radioactive acid waste container. The resin is placed into a resin collection beaker for disposal.

7. Source Materials

7.1. "Isolation of Plutonium, Americium, and Uranium from Urine, by Precipitation, Ion Exchange, and Extraction," 4-16200-RHL-0037, Rocky Flats Plant.

UNITED STATES TRANSURANIUM AND URANIUM REGISTRIES
ANALYTICAL PROCEDURE MANUAL




Figure 1. Ion Exchange Column